

*Surelock*  
**BITAMED**

**BONE MARROW**



*Biopsy System*




# Surelock **BITAMED** BONE MARROW Biopsy System

- Bone marrow punctures are necessary for the diagnosis of various hematological and oncological diseases. As the puncture involves considerable strain on and the risk of infection for the patient.
- The first puncture should ideally yield a biopsy specimen of a sufficient quantity and quality.

**BITAMED** biopsy bone marrow system SURELOCK is fully disposable reducing the risk for infection




- Easy handling , high quality of the materials and workmanship ensure its reliability
- **BITAMED** biopsy bone marrow system SURELOCK features an ergonomic handle allowing for a robust and comfortable grip.
- **BITAMED** biopsy systems were specially designed to facilitate the bone marrow puncture for doctors and to reduce the strain on patients to an absolute minimum.



the Outer needle depth mark allows for accurate targeting and positioning.

Surelock  
**BITAMED** BONE MARROW  
*Biopsy System*




**MARKER :**

► Ensures the size before retrieving the **BITAMED** bone marrow biopsy system SURELOCK. (Every one centimeter . total:4)

**SURELOCK :**

► SURELOCK captures and retrieves quality samples easily and accurately with reduced patient trauma.



**SAFETY CAP :**

► Safety cap also helps the user safety by covering tip after operation.




**Product Ordering Information**

product name	product code	gauge (G)	length (mm)	packing
Bone Marrow biopsy system Surelock	BMB-01	11	100	50 PCS PER BOX Ethylene Oxide Sterile
	BMB-02	13	100	
	BMB-03	11	150	



**BITAMED** Bone marrow biopsy system SURELOCK is also designed for bone marrow aspiration. Connect a luer-lock syringe with 6% luer-taper to the outer needle for bone marrow aspiration.





*Surelock*  
**BITAMED** BONE MARROW  
**Biopsy System**

Factory: NO. 5, Telgerdi 4, Amiriye 3, Ayatollah Rafsanjani highway, Mashhad, IRAN

TeleFax: +9851-35147067

URL: [www.bitamed.ir](http://www.bitamed.ir)

Email: [info@bitamed.ir](mailto:info@bitamed.ir)

Telegram: @bitamed

WhatsApp: +989120436529

Instagram: [daya.daneshgostar](https://www.instagram.com/daya.daneshgostar)



*Scan the barcode to download the digital catalog*



*The death of every green tree is our death*  
**Less paper, more trees**